

# ScopeWriter Instructions

A. Brynes

May 10, 2017

## 1 Introduction

In order to provide an easy method of writing data from the VELA/CLARA LeCroy oscilloscopes to EPICS, a Python program – ScopeWriter – has been written. This document will provide some basic instructions on the functionality of the program. For more information on how the program works, please refer to [1].

## 2 Scope Writer Program

### 2.1 Dependencies

In order to run the Scope Writer program on any LeCroy oscilloscope, the following programs/libraries will be needed:

- Python 2.7.x - <https://www.python.org/>
- PyQt4 - <https://pypi.python.org/pypi/PyQt4>
- PyEPICS - <http://cars9.uchicago.edu/software/python/pyepics3/>
- VELA/CLARA Scope controller - `\\\\fed.cclrc.ac.uk\\\\Org\\\\NLab\\\\ASTeC\\\\Projects\\\\VELA\\\\Software\\\\VELA_CLARA_PYDs\\\\bin\\\\Release`
- PyWin32 - <https://sourceforge.net/projects/pywin32/>

The simplest way to include the Hardware Controller Python libraries to your path is to create an environment variable, `%PYTHONPATH%`, and set it to the network location given above for the VELA/CLARA controllers.

### 2.2 Instructions

- Run `scope_writer_launcher.py`.
- Select controller type ('Physical' for a real scope) and beamline, and click 'Launch'. **NOTE:** At the time of writing (10/5/17), EPICS PVs only exist for the VELA injector, and so this must be selected at the launcher.
- If you wish to use a particular scope setup, load in the setup file using the file browser on the left, and click 'Load setup'. **NOTE:** You will still need to click 'Recall now!' on the oscilloscope interface, under 'File → Recall setup'.
- For each channel you wish to monitor (up to 4 channels can be logged to EPICS at once, but this may slow the loop down to below 10 Hz):

- Select the EPICS PV name to write the calculated value to (these are based on the beamline selected in the launcher).
  - Choose the measurement type (integrated area, maximum/minimum voltage, peak-to-peak).
  - Choose a filter type and interval (currently, only the moving average filter is implemented).
  - Select the region of the scope containing the signal you wish to monitor. **NOTE:** This will depend on the scope setup; currently, the maximum array size we can write to the scope trace PVs is 2000, and so it should be ensured that the timebase of the zoomed-in trace is 1/10 that of the timebase of the main signals. The region of interest that is selected is then a portion of this 2000.
- Click ‘Start logging to EPICS’.
  - To stop the loop, click ‘Abort’.




Figure 1: Launcher for VELA/CLARA ChargeLogger.




Figure 2: GUI for VELA/CLARA ChargeLogger.

## 2.3 Flowchart




Figure 3: Flowchart presenting how the Scope Writer program will work.

## References

- [1] A. D. Brynes. "Trello card: Software: Charge measurement", 2017. <https://trello.com/c/2VQcQ42D/110-software-charge-measurement-adb>.